

**Preventive Maintenance** 

HOSPITAL CONTROL NO.

HOSPITAL NAME			DEPARTMENT			
INSTRUMENT	AutoPulse	MANUFACTURE	Zoll	MODEL		
SERIAL NUMBER						
DATE		Next Due			Period of PM	Months
Action				Pass	Fail	Remark
Visual Inspection th	ne Platform.					
Inspect the Load Plate Cover.						
Inspect the bumpers along the outer edge of the AutoPulse Platform.						
Inspect the head restraint cables.						
Inspect the patient restraint pins.						
Turn the Platform over and inspect the back covers.						
Inspect the two vents on the back of the AutoPulse.						
Inspect the driveshaft.						
Inspect the Battery.	-					
Check the LCD.						
Check the Button.						
The Operation Che	ck					80 PPM ± 5
Cleaning						
Overall Test Resul Comments :	it: PASS / FAIL					

Tested By:
(Signature)

(Name)

**Customer Service**