

448, 450 Soi Supharaj, Phaholyothin Rd., Samsaen Nai, Phyathai, Bangkok 10400, Thailand Tel: 02 090-2591/4 Ext.229 [Service Dept.] Fax: 02 270-0513 [Service Dept.] เลขประจำตัวผู้เลือกาษี: 0105529003833

Preventive Maintenance

HOSPITAL CONTROL NO.

HOSPITAL NAME	DEPARTMENT					
INSTRUMENT	Infusion Pump	MANUFACTURE	Medcaptain	MODEL		
SERIAL NUMBER						
DATE		Next Due			Period of PM	Months
		Action		Pass	Fail	Remark
Overall condition	on / Cleaning					
Mounting						
Power Cord						
Switch, Key Par	nel, LED					
Alarm						
Battery						
Display and Tou	uch screen					
Date & Time Se						
Occlusion Alarm						
Air Bubbler in line						
Pump door Assembly						
Accessory	·					
Connection to S	Station					
Casters						
Electrical Safety						
	•					
Overall Test Resu Comments :	ilt: PASS / FAIL					
				Tested By :		
				(Signature)		
				(Name)		

FP-SEV-05-02 REV.01

Customer Service