

Preventive Maintenance

HOSPITAL CONTROL NO.

HOSPITAL NAME		DEPARTM	MENT		
INSTRUMENT	Patient Monitor	MANUFACTURE Philips	MODEI	1	
SERIAL NUMBER			 -		
DATE		Next Due		Period of PM	Months
			-		
Action			Pass	Fail	Remark
Visual Inspection					
Power On Test					
Basic Performanc	e Assurance Test				
ECG Performance	e Test				
ECG Out Perform	ance Test				
SpO2 Performand	ce Test				
NBP Test					
Invasive Pressure	Performance Test				
Temperature Acc	uracy (in degrees C)				
Mainstream CO2	Accuracy Check				
Sidestream CO2	Accuracy Check				
Cardiac Output Pe	erformance Test				
Power Loss Alarm	n Buzzer Performance Tes	t			
Battery Performar	nce Test				
Touchscreen Cali	bration				
Disabling/Enabling	g Touch Operation				
Print Test Report					
Cleaning					
Overall Test Resul	t: PASS / FAIL				
			Tooled Du	_	
			Tested By (Signature		

FP-SEV-05-02 REV.01

Customer Service

(Name)